**WHO Prequalification of Diagnostics Programme  
PUBLIC REPORT**

**Product: Alere Determine™ HIV-1/2  
Number: PQDx 0033-013-00**

**Abstract**

**Alere Determine™ HIV-1/2** with product codes 7D2342, 7D2343 and 7D2343SET manufactured by **Alere Medical Co. Ltd.**, 357 Matsuhidai Matsudo-shi, Chiba-ken 270-2214, Japan, **rest of the world regulatory version** (non CE-marked regulatory version) was accepted for the WHO list of prequalified diagnostics and was listed on 25 November 2011. This public report was amended on 16 June 2016, and then on 12 July 2016.

**Intended use:**

Alere Determine® HIV-1/2 is an in vitro, visually read, qualitative immunoassay for the detection of antibodies to HIV-1 and HIV-2 in human serum, plasma or whole blood. The test is intended as an aid to detect antibodies to HIV-1/HIV-2 from infected individuals.

**Assay principle:**

Alere Determine HIV-1/2 is an immunochromatographic test for the qualitative detection of antibodies to HIV-1 and HIV-2. Sample is added to the sample pad. As the sample migrates through the conjugate pad, it reconstitutes and mixes with the selenium colloid-antigen conjugate. This mixture continues to migrate through the solid phase to the immobilized recombinant antigens and synthetic peptides at the patient window site.

If antibodies to HIV-1 and/or HIV-2 are present in the sample, the antibodies bind to the antigen-selenium colloid and to the antigen at the patient window, forming a red line at the patient window site. If antibodies to HIV-1 and/or HIV-2 are absent, the antigen-selenium colloid flows past the patient window, and no red line is formed at the patient window site. To ensure assay validity, a procedural control bar is incorporated in the assay device.

A negative result with Determine HIV-1/2 does not exclude the possibility of infection with HIV. A false negative result can occur in the following circumstances:

- low levels of antibody (e.g., early seroconversion specimens) are below the detection limit of the test
- infection with a variant of the virus that is less detectable by the Determine HIV-1/2 assay configuration
- HIV antibodies in the patient that do not react with specific antigens utilized in the assay configuration
- specimen handling conditions which result in loss of HIV antibody multivalency

- For these reasons care should be taken in interpreting negative results. Other clinical data (e.g., symptoms or risk factors) should be used in conjunction with the test results.

Positive specimens should be retested using another method and the results should be evaluated in light of the overall clinical evaluation before a diagnosis is made.

**Test kit contents:**

Component	20 tests (7D2342)	100 tests (7D2343)	100 tests (7D2343SET)
Test cards	2 10-test cards	10 10-test cards	10 10-test cards
Chase buffer	N/A	N/A	1 bottle/2.5ml
EDTA capillary tubes	N/A	N/A	100
Blood lancets	N/A	N/A	100

**Items required but not provided within test kit:**

Consumables	Product codes
Chase buffer (1 bottle/2.5ml)	7D2243
EDTA capillary tubes (100 units)	7D2222
Blood lancets (100 units)	7D2233

**Storage:**

The test kit should be stored at 2-30 °C.

**Shelf-life upon manufacture:**

18 months.

**Summary of prequalification status for Alere Determine™ HIV-1/2**

	Initial acceptance	
	Date	Outcome
<b>PQ amended</b>	16 June 2016, 12 July 2016	listed
<b>Status on PQ list</b>	25 November 2011	listed
<b>Dossier assessment</b>	22 November 2011	MR
<b>Inspection status</b>	24 October 2011	MR
<b>Laboratory evaluation</b>	11 November 2011	MR

MR: Meets Requirements

NA: Not Applicable

**Prioritization for prequalification**

Based on the established criteria, Alere Determine™ HIV-1/2 was given priority for WHO prequalification.

**Product dossier assessment**

Alere Medical Co. Ltd. submitted a product dossier for Alere Determine™ HIV-1/2 as per the “Instructions for compilation of a product dossier” (PQDx\_018 v1). The information submitted in the product dossier was reviewed by WHO staff and external experts (assessors) appointed by WHO.

Commitments for prequalification:

1. analytical performance studies
2. clinical performance studies
3. stability studies
4. a new version of the instructions for use.

WHO will follow-up on implementation of these commitments at the next re-inspection.

Based on the product dossier screening and assessment findings, the product dossier for Alere Determine™ HIV-1/2 meets WHO prequalification requirements.

**Manufacturing site inspection**

A comprehensive inspection was performed at the site of manufacture of Alere Determine™ HIV-1/2 manufactured by Alere Medical Co. Ltd., at 357 Matsuhidai Matsudo-shi, Chiba-ken 270-2214, Japan on 11 to 15 May 2015<sup>1</sup>. The inspection procedure is described in “Information for manufacturers on WHO prequalification inspection procedures for the sites of manufacture of diagnostics” (PQDx\_014 v1).

The inspection found that Alere Medical Co. Ltd. had an established quality management system and manufacturing practices in place that should ensure the manufacture of a product of consistent quality. The manufacturer's responses to the nonconformities found at the time of the inspection were accepted on 13 April 2016.

Commitments for prequalification:

1. Alere Medical Co. Ltd. will continue to review Risk Analysis and Risk Management for accuracy of assessment of risk, attributed to specific components of the product, and the mitigation of such risk, and to ensure ongoing due consideration of end users in resource limited and environmentally challenging regions to which the product is distributed.
2. Alere Medical Co. Ltd. will inform the WHO of changes made subsequent to the site inspection, such as change in location of site of manufacture of major components

---

<sup>1</sup> Previous inspection took place on 28 September to 1 October 2010

of the test, or other changes to the manufacturing process that may affect the quality of the product.

Based on the site inspection and corrective action plan review, the quality management system for Alere Determine™ HIV-1/2 meets WHO prequalification requirements.

### **Laboratory evaluation**

Alere Determine™ HIV-1/2 was evaluated by WHO in the third quarter of 2011 at the Institute of Tropical Medicine, Belgium – a WHO Collaborating Centre for HIV/AIDS Diagnostics and Laboratory Support. The laboratory evaluation was conducted according to the “WHO Protocol for the laboratory evaluation of HIV serology assays” (PQDx\_030 V1.0), and drew the following conclusions:

Alere Determine™ HIV-1/2 is an immunochromatographic rapid diagnostic test for the detection of antibodies to HIV-1/2 in human serum, plasma, and whole blood. A volume of 50µl of serum, plasma or venous/capillary whole blood is required to perform the test procedure. This type of assay requires no sophisticated equipment and can therefore be performed in laboratories with limited facilities. Reading of the results are performed visually i.e. subjective reading.

In this limited performance evaluation using a panel of 1079 biological specimens, we observed an initial (sensitivity (95% CI) of 100% (99.1% - 100%) and an initial specificity (95% CI) of 97.87% (96.4% - 98.8%) compared to the reference assays. The final sensitivity (95% CI) was 100% (99.1% - 100%) and the final specificity (95% CI) was 98.93% (97.8% - 99.6%) compared to the reference assays. In this study, 0.3% of the overall results were recorded as indeterminate. Results were interpreted independently by three technicians; the inter-reader variability was 1.4%. The invalid rate was 0.3% for initial testing and 0.1% for repeat testing.

### **Change notification**

In 2015, Alere Medical Co. Ltd., notified WHO of a change related to addition of a product code for a configuration that included all accessories and reagents required to conduct the test procedure. This notification of change was assessed and product was found to meet WHO prequalification requirements.

In 2016, Alere Medical Co. Ltd., notified WHO of a change related to shelf life. This notification of change was assessed and product was found to meet WHO prequalification requirements.

## **Labelling**

- 1. Labels**
- 2. Instructions for use**

# 1. Labels

Determine™  
**HIV-1/2**

---

**EN**  
For *In Vitro* Diagnostic Use.  
**Alere Determine™ HIV-1/2** is a visual read, qualitative immunoassay for the detection of antibodies to HIV-1 and HIV-2.  
7D2243 Chase buffer is required for whole blood testing.

**Kit contains:**  
• 2 HIV-1/2 recombinant antigen and synthetic peptide coated test cards.

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The Alere Logo and Alere are trademarks of the Alere group of companies.  
© 2015 Alere. All rights reserved.

241106/R6


Determine™  
**HIV-1/2**

---

**EN**  
For *In Vitro* Diagnostic Use.  
**Alere Determine™ HIV-1/2** is a visual read, qualitative immunoassay for the detection of antibodies to HIV-1 and HIV-2.  
7D2243 Chase buffer is required for whole blood testing.

**Kit contains:**  
• 10 HIV-1/2 recombinant antigen and synthetic peptide coated test cards.

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The Alere Logo and Alere are trademarks of the Alere group of companies.  
© 2015 Alere. All rights reserved.

241107/R6

Page 6 of 15

# EDTA Capillary Tubes

List No. 7D2222

100 Each  
*In Vitro*Test

**EN** EDTA Capillary Tubes are to be used in the fingerstick collection of whole blood specimens for use in Determine™ format products. Refer to the assay-specific package insert for full procedure. Bring to room temperature before opening the cap.

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**CAUTION: Glass capillaries may be damaged during transportation or when in use. Handle with care in order to avoid injury when removing from the packaging as well as during use and during disposal. Store in a cool, dry location. Do not refrigerate.**

**Advice Line Telephone number**

Europe & Middle East: + (44) 161 483 9032  
Asia Pacific: + (61) 7 3363 7711  
Africa, Russia & CIS: + (972) 8 9429 683  
Latin America: + (57) 2 6618 797

Store at 2-30°C

©2015 Alere. All rights reserved.  
241203/R6

Produced by:  
Drummond Scientific Co.  
Broomall, PA USA 19008  
For: **Alere Medical Co., Ltd.**

Exp.

Lot

**Alere™ CHASE BUFFER** IVD REF 7D2243  
2.5 mL  
1175831806

LOT XXXXXXXXXXXX  
YYYY-MM-DD  
L060179-01

**Alere Medical Co., Ltd.**  
357 Matsuhidai, Matsudo-shi, Chiba,  
270-2214, Japan

2°C -30°C CE

1135835506

**Alere™ Determine™ Chase Buffer**

**CHASE BUFFER**  
REF 7D2243  
IVD  
Σ 100  
i  
CE 2°C -30°C

**Alere Medical Co., Ltd.**  
357 Matsuhidai,  
Matsudo-shi, Chiba,  
270-2214, Japan  
Tel +81-47-311-5750  
Fax +81-47-311-5751

**EC REP** Alere Ltd.  
Pepper Road,  
Hazel Grove,  
Stockport,  
SK7 5BW, UK

© 2015 Alere. All rights reserved.

135835506

LOT XXXXXXXXXXXX  
YYYY-MM-DD  
L030781-02  
Made in China

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**EN**  
1 Bottle (2,5 mL)  
The Chase Buffer is for use with Determine™ format products.

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**EN**

**Key to symbols used**

<b>REF</b> Catalogue Number	20 Contains Sufficient for 20 tests	Keep away from sunlight
<b>IVD</b> In Vitro Diagnostic Medical Device	100 Contains Sufficient for 100 tests	Do not reuse
Store at 2-30°C	Do not use if package is damaged	

This package insert must be read carefully prior to use. Package insert instructions must be followed accordingly.  
Reliability of assay results cannot be guaranteed if there are deviations from the instructions in this package insert.

**NAME AND INTENDED USE**  
The **Alere Determine™ HIV-1/2** is an *In Vitro*, visually read, qualitative immunoassay for the detection of antibodies to HIV-1 and HIV-2 in human serum, plasma or whole blood. The test is intended as an aid to detect antibodies to HIV-1/HIV-2 from infected individuals.

**SUMMARY AND EXPLANATION OF THE TEST**  
AIDS (Acquired Immunodeficiency Syndrome) is characterized by changes in the population of T-cell lymphocytes. In an infected individual, the virus causes depletion of helper T-cells, which leaves the person susceptible to opportunistic infections and some malignancies. The virus that causes AIDS exists as two related types known as HIV-1 and HIV-2. The presence of the AIDS virus elicits the production of specific antibodies to either HIV-1 or HIV-2.<sup>1,2</sup>

**BIOLOGICAL PRINCIPLES OF THE PROCEDURE**  
**Alere Determine™ HIV-1/2** is an immunochromatographic test for the qualitative detection of antibodies to HIV-1 and HIV-2.  
Sample is added to the sample pad. As the sample migrates through the conjugate pad, it reconstitutes and mixes with the selenium colloid-antigen conjugate. This mixture continues to migrate through the solid phase to the immobilized recombinant antigens and synthetic peptides at the patient window site. If antibodies to HIV-1 and/or HIV-2 are present in the sample, the antibodies bind to the antigen-selenium colloid and to the antigen at the patient window, forming a red line at the patient window site. If antibodies to HIV-1 and/or HIV-2 are absent, the antigen-selenium colloid flows past the patient window, and no red line is formed at the patient window site.  
To insure assay validity, a procedural control bar is incorporated in the assay device.

**CONTENTS**  
**Alere Determine™ HIV-1/2 Serum/Plasma Assay, 20 Tests (7D2342) or 100 Tests (7D2343)**

- **Alere Determine™ HIV-1/2 Test Card**, 2 cards or 10 cards (10 tests/card), HIV-1/2 recombinant antigen and synthetic peptide coated.

**ACCESSORIES (required but not provided)**  
**For testing Whole Blood samples**

- 1 Bottle (2.5 mL) Chase Buffer (7D2243) prepared in phosphate buffer. Preservatives: Antimicrobial Agents.

**Whole Blood (Fingerstick assay)**

- EDTA Capillary Tubes (7D2222)

**Materials Required But Not Provided**  
Disposable gloves, timing device  
Micropipette capable of delivering 50 µL (other than fingerstick)  
Alcohol swab, gauze pad, Lancet (for fingerstick)

**WARNINGS AND PRECAUTIONS**  
**For In Vitro Diagnostic Use.**

**CAUTION:**  
Appropriate biosafety practices<sup>3</sup> should be used when handling specimens and reagents. These precautions include, but are not limited to the following:

- Wear gloves.
- Do not pipette by mouth.
- Do not eat, drink, smoke, apply cosmetics, or handle contact lenses in areas where these materials are handled.
- Clean and disinfect all spills of specimens or reagents using a suitable disinfectant, such as 0.5% sodium hypochlorite.<sup>4,7</sup>
- Decontaminate and dispose of all specimens, reagents, and other potentially contaminated materials in accordance with local regulations.<sup>4,8</sup>

**STORAGE**  
The **Alere Determine™ HIV-1/2 Test Cards** and **Chase Buffer** must be stored at 2-30°C until expiration date. Kit components are stable until expiration date when handled and stored as directed. Do not use kit components beyond expiration date.  
Immediately reseal all unused tests in the foil pouch containing the desiccant by pressing seal from end to end to close.

**SPECIMEN COLLECTION**  
**Serum, Plasma, and Whole Blood Collection by Venipuncture**  
Human serum, plasma, and whole blood collected by venipuncture should be collected aseptically in such a way as to avoid hemolysis.

**NOTE: For whole blood and plasma specimens, EDTA collection tubes must be used.**

**Whole Blood Collection by Fingerstick<sup>9</sup>**  
Before collecting a fingerstick specimen, place an EDTA capillary tube on a clean dry surface.

1. Choose the fingertip of the middle, ring, or index finger (whichever is the least callused) for adults and children older than one year. Warm the hand as needed with a warm, moist towel or warm water to increase blood flow.
2. Clean fingertip with alcohol; allow to air dry. Position the hand palm-side up.
3. Use a new lancet for each person. Place the lancet off-center on the fingertip. Firmly press the lancet against the finger and puncture the skin. Dispose of the lancet in an appropriate biohazard sharps container.
4. Wipe away the first drop of blood with a sterile gauze pad.
5. Hold the finger lower than the elbow and apply gentle, intermittent pressure to the base of the punctured finger several times. Touch the tip of the EDTA Capillary Tube to the drop of blood\*. Avoid air bubbles.

\*If EDTA Capillary Tubes (7D2222) will be used, fill the tube with blood between the 2 marked lines (50 µL).

**SPECIMEN STORAGE**

- Serum and plasma specimens should be stored at 2-8°C if the test is to be run within 7 days of collection. If testing is delayed more than 7 days, the specimen should be frozen (-20°C or colder).
- Whole blood collected by venipuncture should be stored at 2-8°C if the test is to be run within 7 days of collection.  
Do not freeze whole blood specimens.
- Whole blood collected by fingerstick should be tested immediately.

**TEST PROCEDURE**  
The desired number of test units from the 10-test card can be removed by bending and tearing at the perforation.

**NOTES:**

- Removal of the test units should start from the right side of the test card to preserve the lot number which appears on the left side of the test card.
- Assay should be initiated within 2 hours after removing the protective foil cover from each test.

1. Remove the protective foil cover from each test.
2. For serum or plasma samples:
  - a. Apply 50 µL of sample (precision pipette) to the sample pad (marked by the arrow symbol).
  - b. Wait a minimum of 15 minutes (up to 60 minutes) and read result.
3. For whole blood (venipuncture) samples:
  - a. Apply 50 µL of sample (precision pipette) to the sample pad (marked by the arrow symbol).
  - b. Wait one minute, then apply one drop of Chase Buffer to the sample pad.
  - c. Wait a minimum of 15 minutes (up to 60 minutes) and read result.
4. For whole blood (fingerstick) samples:
  - a. Apply 50 µL of sample (by EDTA capillary tube) to the sample pad (marked by the arrow symbol).
  - b. Wait until blood is absorbed into the sample pad, then apply one drop of Chase Buffer to the sample pad.
  - c. Wait a minimum of 15 minutes (up to 60 minutes) and read result.

**QUALITY CONTROL**

To insure assay validity, a procedural control is incorporated in the device and is labeled "Control". If the control bar does not turn red by assay completion, the test result is invalid and the sample should be retested.

**INTERPRETATION OF RESULTS**

**POSITIVE (Two Bars)**

Red bars appear in both the control window (labeled "Control") and the patient window (labeled "Patient") of the strip. Any visible red bar in the patient window should be interpreted as positive.


**NEGATIVE (One Bar)**

One red bar appears in the control window of the strip (labeled "Control"), and no red bar appears in the patient window of the strip (labeled "Patient").

**INVALID (No Bar)**

If there is no red bar in the control window of the strip, and even if a red bar appears in the patient window of the strip, the result is invalid and should be repeated.

**NOTES:**

- The test result is positive even if the patient bar appears lighter or darker than the control bar.
- The control bar may exhibit a weak intensity for some patient samples, particularly those with high titer HIV.
- Upon appearance of a red bar in the control window, no matter how faint, the test result is considered valid.
- If an invalid test result occurs repeatedly, or for technical assistance, contact your local distributor or call Technical Support.

**LIMITATIONS OF THE PROCEDURE**

- The **Alere Determine™ HIV-1/2** test is designed to detect antibodies to HIV-1 and HIV-2 in human serum, plasma, and whole blood. Other body fluids or pooled specimens may not give accurate results.
- The intensity of the patient bar does not necessarily correlate to the titer of antibody in the specimen.
- A negative result with **Alere Determine™ HIV-1/2** does not exclude the possibility of infection with HIV. A false negative result can occur in the following circumstances:
  - low levels of antibody (e.g., early seroconversion specimens) are below the detection limit of the test
  - infection with a variant of the virus that is less detectable by the **Alere Determine™ HIV-1/2** assay configuration
  - HIV antibodies in the patient that do not react with specific antigens utilized in the assay configuration
  - specimen handling conditions which result in loss of HIV antibody multivalency
  - HIV-infected persons taking antiretroviral medication<sup>11,12,13</sup>

For these reasons care should be taken in interpreting negative results. Other clinical data (e.g., symptoms or risk factors) should be used in conjunction with the test results.

- Positive specimens should be retested using another method and the results should be evaluated in light of the overall clinical evaluation before a diagnosis is made.
- Whole blood or plasma specimens containing anticoagulants other than EDTA may give incorrect results.
- Neonates of HIV-infected mothers may carry maternal antibodies to HIV for up to around eighteen months, which may not necessarily indicate the true infection status of the new born.

**PERFORMANCE CHARACTERISTICS**

**SPECIFICITY**

A total of 1594 serum and plasma specimens from Asia, West Africa, and North America were tested by **Alere Determine™ HIV-1/2** and a commercially available test (Table I).

**Table I**  
Specificity of **Alere Determine™ HIV-1/2**

Population	Number of Specimens Tested	Negative by <b>Alere Determine™ HIV-1/2</b>	Negative by a Commercially Available Test***
Seronegative Serum	908	907/908 (99.89%)	908/908 (100.00%)
Plasma	403	403/403 (100.00%)	403/403 (100.00%)
Pregnant Females	58*	57/57 (100.00%)	57/57 (100.00%)
West Africans	49	48/49 (97.96%)	48/49 (97.96%)
Disease States Other than HIV and Potentially Interfering Substances	176*	173/175 (98.86%)	174/175 (99.45%)
<b>Total</b>	<b>1594**</b>	<b>1588/1592 (99.75%)</b>	<b>1590/1592 (99.87%)</b>

\* One specimen from a pregnant female and an HCV positive patient were positive by both **Alere Determine™ HIV-1/2** and the commercially available test. Both specimens confirmed positive by HIV-1 Western Blot.

\*\* 456 specimens were from North America, 1089 specimens were from Asia, and 49 specimens were from Africa.

\*\*\* The reference method of a commercially available test is particle agglutination.

A total of 3663 seronegative serum and plasma specimens from North America, Asia, and Africa were tested by **Alere Determine™ HIV-1/2** and commercially available tests (Table II). The specimens from North America, Asia, and 49 of 2118 specimens from Africa (referred to as "West Africans" in Table I) were included in Table I. Discordant specimens were confirmed negative by either Western blot or HIV-1 PCR assays.

**Table II**  
A Comparison of **Alere Determine™ HIV-1/2** Specificity by Geographic Area

Area	Number of Specimens Tested	Negative by <b>Alere Determine™ HIV-1/2</b>	Negative by Commercially Available Tests*
North America	456	451/454 (99.34%)	453/454 (99.78%)
Asia	1089	1089/1089 (100.00%)	1089/1089 (100.00%)
Africa	2118	2079/2118 (98.16%)	2100/2118 (99.15%)

\* The reference methods of commercially available tests are particle agglutination, enzyme immunoassay and chemiluminescent immunoassay.

A total of 368 seronegative whole blood specimens from Thailand were tested with paired serum and plasma by **Alere Determine™ HIV-1/2**. Thirty-nine of the whole blood specimens were collected by both venipuncture and fingerstick (Table III).

**Table III**  
A Comparison of **Alere Determine™ HIV-1/2** Specificity in Seronegative Whole Blood and Paired Serum and Plasma Specimens

Specimen Type	Number of Specimens Tested	Negative by <b>Alere Determine™ HIV-1/2</b>
Serum	368	368/368 (100.00%)
Plasma	368	368/368 (100.00%)
Whole Blood (venipuncture)	368	368/368 (100.00%)
Whole Blood (fingerstick)	39	39/39 (100.00%)

**SENSITIVITY**

A total of 869 HIV-1 and HIV-2 antibody positive serum and plasma specimens from Asia, Africa, North and South America were tested by **Alere Determine™ HIV-1/2** and a commercially available test (Table IV).

**Table IV**  
Sensitivity of **Alere Determine™ HIV-1/2**

Population	Number of Specimens Tested	Positive by <b>Alere Determine™ HIV-1/2</b>	Positive by a Commercially Available Test**
HIV-1 Positive	521*	521/521 (100.00%)	521/521 (100.00%)
HIV-2 Positive	114*	114/114 (100.00%)	114/114 (100.00%)
HIV-1 Subtypes A-G	222	222/222 (100.00%)	Not Tested Not Tested
HIV-1 Group O	12	12/12 (100.00%)	Not Tested Not Tested
<b>Total</b>	<b>869</b>	<b>869/869 (100.00%)</b>	<b>635/635 (100.00%)</b>

\* 228 specimens were from North America, 296 specimens were from Asia, and 111 specimens were from Africa.

\*\* The reference method of a commercially available test is particle agglutination.

A total of 1653 seropositive serum and plasma specimens from North America, Asia, and Africa were tested by **Alere Determine™ HIV-1/2** and commercially available tests (Table V). The specimens from North America, Asia, and 111 of 1129 specimens from Africa (referred to as "HIV-2 Positive" in Table IV) were included in Table IV. Discordant specimens were confirmed HIV-1 positive by either Western blot or HIV-1 PCR assays.

**Table V**  
A Comparison of **Alere Determine™ HIV-1/2** Sensitivity by Geographic Area

Area	Number of Specimens Tested	Positive by <b>Alere Determine™ HIV-1/2</b>	Positive by Commercially Available Tests**
North America	228	228/228 (100.00%)	228/228 (100.00%)
Asia	296	296/296 (100.00%)	296/296 (100.00%)
Africa	1129	1128/1129 (99.81%)	1129/1129 (100.00%)

\* One negative specimen by **Alere Determine™ HIV-1/2** confirmed positive by HIV-1 PCR.

\*\* The reference methods of commercially available tests are particle agglutination, enzyme immunoassay and chemiluminescent immunoassay.

A total of 102 seropositive whole blood specimens from Thailand were tested with paired serum and plasma by **Alere Determine™ HIV-1/2**. Thirty-two of the whole blood specimens were collected by both venipuncture and fingerstick (Table VI).

**Table VI**  
A Comparison of **Alere Determine™ HIV-1/2** Sensitivity in Seropositive Whole Blood and Paired Serum and Plasma Specimens



Specimen Type	Number of Specimens Tested	Positive by <b>Alere Determine™ HIV-1/2</b>
Serum	102	102/102 (100.00%)
Plasma	102	102/102 (100.00%)
Whole Blood (venipuncture)	102	102/102 (100.00%)
Whole Blood (fingerstick)	32	32/32 (100.00%)

The manufacturing process produces different lot numbers for the kit and test cards; these lot numbers are traceable.


**BIBLIOGRAPHY**

- Plot P, Plummer FA, Mhalu FS, Lamborany JL, Chin J, Mann JM. AIDS: An International Perspective. *Science*. 1988; 239: 573-579.
- Weniger BG, Takebe Y, Ou CY, et al. The Molecular Epidemiology of HIV in Asia. *AIDS*. 1994; 8(S2): S13-S28.
- Gürtler LG, Hauser PH, Eberle J, et al. A New Subtype of Human Immunodeficiency Virus Type 1 (MVP-5180) from Cameroon. *Journal of Virology*. 1994; 68(3): 1581-1585.
- World Health Organization. Laboratory Biosafety Manual. Geneva: World Health Organization, 1993.
- National Committee for Clinical Laboratory Standards. *Protection of Laboratory Workers from Infectious Disease Transmitted by Blood, Body Fluids, and Tissue*. Tentative Guideline. NCCLS Document M29-12. Villanova, PA: NCCLS, 1991: 1-43.
- CDC. Recommendations for Prevention of HIV Transmission in Health-Care Settings. *MMWR* 1987; 36(2S): 3S-18S.
- Sehuster LM, Hollinger FB, Dresman GR, et al. Immunological and Biophysical Alteration of Hepatitis B Virus Antigens by Sodium Hypochlorite Disinfection. *Applied and Environmental Microbiology*. 1981; 42(5): 762-7.
- Clinical and Laboratory Standards Institute. *Clinical Laboratory Waste Management: Approved Guideline - Third Edition*. GP05-A3 Vol.31 No.3 January 2011
- EPA Guide for Infections Waste Management: Publication No. EPA/530-SW-86-014. Washington, DC: US Environmental Protection Agency, 1986:1-1 - 5.5, R1-R3, A1-A24.
- Clinical and Laboratory Standards Institute. *Procedures and Devices for the Collection of Diagnostic Capillary Blood Specimens*; Approved Standard-Sixth Edition GP42-A6 Vol.28 No.25 September 2008
- Delaney KR, Branson BM, et al. Evaluation of the Performance Characteristics of 6 Rapid HIV Antibody Tests. *Clinical Infectious Diseases*. 2011; 52(2): 257-263.
- O'Connell RJ, Merritt TM, Malia JA, et al. Performance of the OraQuick Rapid Antibody Test for Diagnosis of Human Immunodeficiency Virus Type 1 Infection in Patients with Various Levels of Exposure to Highly Active Antiretroviral Therapy. *Journal of Clinical Microbiology*. 2003; 41(5):2153-2155.
- O'Connell RJ, Ajan BK, Anderson SA, Malia JA, Michael NL. Sensitivity of the Multispot HIV-1/HIV-2 Rapid Test Using Samples from Human Immunodeficiency Virus Type 1-Positive Individuals with Various Levels of Exposure to Highly Active Antiretroviral Therapy. *Journal of Clinical Microbiology*. 2006; 44(5): 1831-1833.

**Whole Blood**

**Serum, Plasma**

**Advice Line**

For further information, please contact your distributor, or call to one of the following Alere Product Support Care Centers:

Region	Phone	E-Mail Address
Europe & Middle East	+ (44) 161 483 9032	EMproductsupport@alere.com
Asia Pacific	+ (61) 7 3363 7711	APproductsupport@alere.com
Africa, Russia & CIS	+ (972) 8 9429 683	ARCISproductsupport@alere.com
Latin America	+ (57) 2 661 8797	Lproductsupport@alere.com

**Alere Medical Co., Ltd.**  
 357 Matsuhidai, Matsudo-shi,  
 Chiba, 270-2214 Japan  
[www.alere.com](http://www.alere.com) [www.alerehiv.com](http://www.alerehiv.com) Tel +81 47 311 5750


The Alere Logo, Alere and Determine are trademarks of the Alere group of companies.  
 © 2015 Alere. All rights reserved.

## Determine™ Chase Buffer

REF 7D2243

Key to symbols used/Používané symboly/Symbolforklaring/Erläuterung der verwendeten Symbole/  
Πίνακας συμβόλων/Clave de los símbolos utilizados/Légende des symboles utilisés/  
Legenda dei simboli utilizzati/Legenda dos símbolos utilizados/Symbolforklaring

Store at 2-30°C/Skladujte při teplotě 2-30°C/  
Opbevar ved 2-30°C/Lagerung bei 2 bis 30°C/  
Φυλάσσεται στους 2-30°C/Almacenar a 2-30°C/  
Conserver entre 2 et 30°C/Conservare a 2-30°C/  
Conservar a 2°C-30°C/Förvaras vid 2-30°C

### CHASE BUFFER

Chase Buffer/Nosný pufr/  
Chase-buffer/Antriebspuffer/  
Ρυθμιστικό διάλυμα σταθεροποίησης/  
Tampón de arrastre/Tampón de fixation/  
Tampone chase/Tampão de fixação/  
Fixeringsbuffer

For In Vitro Diagnostic Use/Pro diagnostické účely in vitro/  
Til in vitro-diagnostisk brug/  
Der Test ist nur für die In-vitro-Diagnos vorgesehen/  
Για in vitro διαγνωστική χρήση/Para uso en diagnósticos in vitro/  
Pour usage diagnostic in vitro/Per uso diagnostico In vitro/  
Para uso em Diagnóstico In Vitro/För diagnostisk användning in vitro

## EN

### Name and Intended Use

The Chase Buffer is for use with Determine™ format products.  
Refer to the assay-specific package insert for additional information.  
When adding Chase Buffer to the sample pad, hold the bottle vertically.  
One bottle of Chase Buffer can be used for 100 tests.

### Contents

**CHASE BUFFER** 1 Bottle (2.5 mL) Chase Buffer prepared in phosphate buffer.  
Preservatives: Antimicrobial Agents.

### Storage Instruction

Recap and store the chase buffer at 2-30°C to avoid evaporation or spillage.

### Advice Line

For further information, please contact your distributor, or call Aleré Technical Specialists:

Africa, Russia & CIS: Tel: +972 8 9429 683  
Email: ARCISproductsupport@alere.com

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